



ELISA Technologies, Inc.

elisa-tek.com



PAGE 1 OF 1 - REPORT 2104117

REPORT DATE: April 22, 2021

REPORT TO: *Dominic Reale
Goosefoot Acres Inc.
5879 Center Rd
Valley City, OH 44280
Tel: 440-476-3802
Email: dreale.dandyblend@gmail.com*

SUBJECT: *PO Number: Dominic; Sample for Gliadin Analysis*

RECEIVED: *April 16, 2021;
USPS Priority Mail 3-Day TRK# 9405 5118 9956 0147 8593 94*

ANALYTICAL RESULTS: RESULTS APPLY ONLY TO THE SAMPLES AS RECEIVED.

CONTROL #	SAMPLE DESCRIPTION	ANALYSIS IDENTIFICATION	Gliadin
2104117-1	<i>Dandy Blend 14oz Lot# 0521A Ex. Date 12/21/2023</i>	500828 <i>Gliadin</i> Assay	Not Detected

NOTE: The RIDASCREEN FAST Gliadin assay is a sandwich ELISA for the quantitative analysis of prolamins from wheat (gliadin), rye (secalin), and barley (hordein) in food. Results are expressed in ppm (mg/kg) gliadin with a quantifiable range of 5.0-40.0 ppm. A result of "not detected" indicates that gliadin was not detected in the sample (<5.0 ppm). To obtain gluten concentrations, it is acceptable to multiply the gliadin value by two.

Information in italics has been supplied by the customer.

Respectfully Submitted for ELISA Technologies, Inc.

By 

Analyst

By 

Eun Park, LFA Director

End of Report

THIS REPORT SHALL NOT BE REPRODUCED EXCEPT IN ITS ENTIRETY AND WITH THE WRITTEN PERMISSION OF ELISA TECHNOLOGIES, INC.